

JH-W2





Product feature

- 1. Modern & ITE design, discreet and light with only 5g;
-Powerful hearing support in a wireless earbud form factor
- 2. 5.0 bluetooth function, faster and more stable connection
- 3. Clear and fidelity sound with DSP chip
-4 Channels /12 bands / WDRC /AFC /NR features
- 4. Calls and Music auto switch
-Binaural Bluetooth connection & Call and stream via Bluetooth
- 5. Fast charge technology
-Portable charging case can recharge double devices at least by 5 times
-Quick 15-minute charge provides 8+hours hearing back.
-20+ Hours of hearing aid Use,3 Hours of Bluetooth streaming
-2 Hours full charge for charging base

Technical Data	IEC 60118-7:2005 Test Data (2cc-coupler)	ANSI S3.22-2003 Test Data (2cc-coupler)
Operating Voltage	DC 3.0V~4.25V	DC 3.0V~4.25V
Acoustic Gain (50 dB SPL)		
HFA	26.3dB	22.9dB
1600Hz	25dB	24dB
Peak Value	28.9dB	/
Max. output		
HFA	108.3dB	107.3dB
1600Hz	104dB	107dB
Peak Value	110.6dB	109.1dB
Reference Test Gain	24dB	25.2dB
Induction Coil Sensitivity	/	/
Frequency Range	216~4292Hz	200-4290Hz
Total Harmonic Distortion(500/800/1600Hz)	0.5%/0.2%/0.3%	0.4%/0.8%/0.2%
Equivalent Input Noise	25.8dB	28dB
Battery Current	1.85mA	1.85mA
Battery Type	Lithium battery	Lithium battery
Working Time	20h*5 with charging base	20h*5 with charging base

FITTING RANGE



- 1. Hearing aid charging indicator /Bluetooth indicator LED light
- 2. Volume Control
- 3. Bluetooth Microphone
- 4. Hearing aid Microphone
- 5. Multi-function button (Bluetooth switch)
- 6. Charging Contacts

JH-W2

Accessories



Ear Dome(L)



Ear Dome(M)



Ear Dome(S)




Micro USB Cable




Cleaning tool

All Curves are measured with fonix 8000 hearing aid test system, 2ccm-coupler ,IEC 60118-7.All sound pressure levels are referred to 20μPa (OSPL90 Output/FOG50 Gain/Reference Test Gain(RTG))


OSPL₉₀ Output



FOG₅₀ Gain



Reference Test Gain (RTG)



WARNING! Choking hazard posed by small parts.
Infants, small children and persons of mental incapacity
must not wear the hearing instrument
without appropriate supervision.