

PRODUCT FEATURES


- * Multiple directionality microphones and digital signal processing (channels)
- * 4 Hearing programs
- * Binaural type
- * Hearing App (iOS and Android)
- * Direct Audio Streaming from iOS devices or Android devices
- * Wireless programming
- * Feedback cancellation & WDRC
- * Multi-function Button
- * Rechargeable




ANSI 3.22 2014 2cc coupler technical data

Test Item	Mode One	Mode Two	Mode Three	Mode Four
Max. Output				
Peak Value	104dB SPL	105dB SPL	104dB SPL	109dB SPL
Gain				
HFA-FOG50	23dB	19dB	25dB	19dB
The self- generated noise value(EIN)	15dB	21dB	10dB	16dB
Frequency Range	210Hz-7300Hz	220Hz-7500Hz	220Hz-6900Hz	200Hz-7500Hz
Total Harmonic Distortion	0.3%	0.3%	0.3%	0.2%
Charging Time	≤1h			
Battery Type	Rechargeable Li-ion Battery			
Working Time	Up to 15 hours (Bluetooth disconnected)			
	Up to 10 hours (Bluetooth connected but not use streaming and App)			
	Up to 2 hours of use with streaming (e.g: phone call and music)			
	The fully charging case supports to 3X charging cycles for a pair of hearing aids.			

***The data sheet report is only responsible for the sample**



PART OF HEARING AID



1. Microphone
2. Indicator Light
3. Multi-function Button
(On/Off, switches hearing programs, adjusts volume level)
4. Speaker Unit
5. Eartip
6. Thin Tubing

WARNING! Choking hazard posed by small parts. Infants, small children and persons of mental incapacity must not wear the hearing instrument without appropriate supervision.

***The data sheet report is only responsible for the sample**