

JH-909G





Product feature

- 1.Amplifies the volume of all sounds indoors and outdoors
- 2.With micro USB charger
- 3.Rechargeable hearing amplifiers
- 4.Ergonomic fit directly in the ear canal;
- 5.Single ear use,suitable for left or right ear

Technical Data	IEC 60118-7:2005 Test Data (2cc-coupler)	ANSI S3.22-2003 Test Data (2cc-coupler)
Operating Voltage	DC 3.0V~4.25V	DC 3.0V~4.25V
Acoustic Gain (50 dB SPL)		
HFA	27.2dB	27.1dB
1600Hz	28.9dB	28.5dB
Peak Value	28.9dB	N/A
Max. output		
HFA	110.1dB	110.2dB
1600Hz	113.2dB	112.9dB
Peak Value	115.3dB	115.4dB
Reference Test Gain	27.1dB	27.1dB
Induction Coil Sensitivity	N/A	N/A
Frequency Range	283Hz-3921Hz	283Hz-3934Hz
Total Harmonic Distortion(500/800/1600Hz)	1%/1.8%/0.3%	1.4%/1.6%/0.3%
Equivalent Input Noise	26.9dB	26.6dB
Battery Current	1.1mA	1.1mA
Battery Type	Lithium battery	Lithium battery
Working Time	16*3 with charging base	16*3 with charging base

FITTING RANGE



1. Microphone
2. Multi-function Button (ON/OFF, Volume Control)
3. Ear Dome
4. Charging Contacts

JH-909G

Accessories



Ear Dome(L)



Ear Dome(M)



Ear Dome(S)




Micro USB Cable




Cleaning tools

All Curves are measured with fonix 8000 hearing aid test system, 2cc-coupler ,IEC 60118-7.All sound pressure levels are referred to 20μPa (OSPL90 Output/FOG50 Gain/Reference Test Gain(RTG))


OSPL₉₀ Output



FOG50 Gain



Reference Test Gain (RTG)



WARNING! Choking hazard posed by small parts.
Infants, small children and persons of mental incapacity
must not wear the hearing instrument
without appropriate supervision.

*The data sheet report is only responsible for the sample