


# JH-338-C





## Product feature

- 1.Rechargeable  
-Micro USB charging dock , 18 hours of hearing on one charge
- 2.Stable and Reliable  
-Classic Design with Good quality
- 3.Easy Operation  
-Simple slider key operation for volume  
-Classic rocker-key operation to turn on and off, more adapt to the operation habits of the elderly
- 4.Moving coil loudspeaker in the ear, great gain;
- 5.Rotated tube, fit for either ear;

| Technical Data                            | EN 60118-7: 2005<br>(2 ccm-coupler) | ANSI S3.22-2003<br>(2 ccm-coupler) |
|-------------------------------------------|-------------------------------------|------------------------------------|
| Operating Voltage                         | DC1.0V~1.5V                         | DC1.0V~1.5V                        |
| Acoustic Gain (50 dB SPL)                 |                                     |                                    |
| HFA                                       | 29.5dB                              | 29.7dB                             |
| 1600Hz                                    | 32.5dB                              | 32dB                               |
| Peak Value                                | 34.8dB                              | N/A                                |
| Max. output                               |                                     |                                    |
| HFA                                       | 111.6dB                             | 111.6dB                            |
| 1600Hz                                    | 113dB                               | 112dB                              |
| Peak Value                                | 116dB                               | 116dB                              |
| Reference Test Gain                       | 29.6dB                              | 29.7dB                             |
| Induction Coil Aensitivity                | /                                   | /                                  |
| Frequency Range                           | 345Hz-3894Hz                        | 344Hz-3890Hz                       |
| Total Harmonic Distortion(500/800/1600Hz) | 1.5%/0.8%/0.2%                      | 2.0%/0.6%/0.1%                     |
| Equivalent Input Noise                    | 28dB                                | 28dB                               |
| Battery Current                           | 0.95mA                              | 0.95mA                             |
| Battery Type                              | NI-MH battery                       | NI-MH battery                      |
| Working Time                              | 18h                                 | 18h                                |




## FITTING RANGE



- 1.Microphone
- 2.Volume Control
- 3.O/N
- 4.Charging Contacts
- 5.Ear Tube
- 6.Ear Dome


# JH-338-C

## Accessories




All Curves are measured with fonix 8000 hearing aid test system, 2ccm-coupler ,IEC 60118-7.All sound pressure levels are referred to 20μPa (OSPL90 Output/FOG50 Gain/Reference Test Gain(RTG))


OSPL<sub>90</sub> Output



FOG<sub>50</sub> Gain



Reference Test Gain (RTG)



**WARNING!** Choking hazard posed by small parts.  
Infants, small children and persons of mental incapacity  
must not wear the hearing instrument  
without appropriate supervision.